## Descripción

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## Instructions for use

- Warnings**
- This product is supplied sterile and are intended for single use only. Do not refill, re-use and re-sterilize.
  - Do not use the Syrinjector for the administration of lipid emulsions or medicinal solutions containing lipid emulsion; medicinal solutions containing any oily ingredients; surface-active agents; highly viscous medicinal solutions; or medicinal solutions containing solubilizing agents such as ethanol.

### Precautions

- Set the flow rate by taking into consideration the target site, patient's physique and age, user's experience and technique.
- Avoid tightening the joint parts excessively with tightening tools, etc. The product may crack, and the drug solution may leak out.
- Follow the established procedure when using the product, and always expel any air that collects in the Syrinjector.
- When filling the product with medicinal solution, make sure it is free of glass shards or other contaminants.
- Operate the product slowly when filling.
- Do not subject the infusion port to excessive strain when filling.
- Fill the product under sterile conditions.
- Fill the medicinal solution to the prescribed capacity (60 mL or 120 mL).
- Cap the infusion port after filling.
- Use the product immediately after filling.
- Prevent additionally administered medicinal fluid from flowing toward the flow selector when using the 3-way stopcock between the flow selector and catheter
- Make sure that organic solvents such as alcohol do not come into contact with the product, or connectors.
- The flow rate is most accurate when saline solution is used at 23°C room temperature. The viscosity of the solution may be affected by the temperature of the solution (drug &/or diluent), the concentration of the solution, the pressure thereby impacting the flow rate.
- Note that flow rate may be affected by factors such as the insertion site and length or inner diameter of connected catheters.
- A trace amount of medicinal solution may remain in the pump or tubes.
- Do not use the product in high-pressure oxygen.
- Using lipophilic drugs presents a risk of elution of bis(2-ethylhexyl) phthalate, a PVC plasticizer.
- If a narcotic drug is used in a product of the type in which the medicinal solution can be removed through the injection port, it must only be handled within a medical facility.
- Do not use the product with emulsive or highly viscous medicinal solutions.
- When the product is used with a nerve block or subcutaneous injection, the amount ejected may be less than the prescribed amount, or injection may stop. Also be aware of back-pressure from arterial injection.
- Before using the product, examine the connection status of catheters, tubes and connectors, check for any tube breakage, etc., and confirm that the product will operate normally. Re-check the catheters to be placed inside the patient's body.
- When using a medicine that may crystallize, check regularly if the injection is performed appropriately, when using a medical material that could be crystallized.

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<b>How to Use</b>	<b>ENG</b>
<ol style="list-style-type: none"><li>Remove the Syrinjector from its sterilized package.</li> <li>i) Make sure that the Flow Selector is adjusted to the maximum flow rate. If not, please use the hexagonal slot lever on the back of the case to adjust flow rate to the maximum value. (A) <ul style="list-style-type: none"><li>ii) Open the 3-way cock to the straight direction only.(B)</li></ul></li> <li>Open the cap of the injection port. (C)</li> <li>Insert the tip of the filled syringe into the injection port and turn to lock. Inject the medicinal solution slowly under aseptic condition. Mixture of glass tip or other substances in a medicinal solution may cause the leak from injection port.</li> <li>After the medicinal solution has been injected, remove the syringe from the Syrinjector by twisting it off.</li> <li>Close the cap of the injection port.</li> <li>Remove the Priming cap, and make sure that the medicinal solution is flowing out of the patient-end of the Flow Selector or the 3-way cock.</li> <li>Set the indicator of Flow Selector to the "OFF" position. (A)</li> <li>Connect the patient-end of the Flow Selector or the 3-way cock to the indwelling catheter.</li> <li>Set the Flow Selector to the target flow rate. (A)</li> <li>Write the necessary information on the Name tag and attach it to the Syrinjector.</li> <li>Put the Syrinjector into the carrying bag if necessary.</li></ol>	

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